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The EDV28A-N-01 is an in-place vital signs monitoring sensor based on 60GHz millimeter-wave radar technology. It can detect millimeter-level chest movements in real time, providing comfortable, non-contact long-term stable monitoring. Combining edge computing and intelligent algorithms, it continuously monitors and analyzes the target's presence status as well as key indicators such as heart rate, respiration, and body movement.

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Non-contact Monitoring: Utilizing cutting-edge 60GHz millimeter-wave sensing technology to detect millimeter-level chest movements in real time, providing comfortable, non-contact long-term stable monitoring.

- 智能分析：结合边缘计算和智能算法，连续监测分析目标在位情况，及心率、呼吸、体动等关键指标。
- Intelligent Analysis: Combining edge computing and intelligent algorithms to continuously monitor and analyze the target's presence status, heart rate, respiration, and body movement.

- 睡眠监测：智能监测睡眠阶段的浅睡、深睡、REM 睡眠、清醒期等情况，可评估和复盘睡眠质量。
- Sleep Monitoring: Intelligently monitors sleep stages including light sleep, deep sleep, REM sleep, and

wakefulness, enabling sleep quality assessment and review.

- 异常警报：智能识别异常呼吸或心率，支持及时警报，守护用户健康安全。
Abnormal Alert: Intelligently identifies abnormal respiration or heart rate, supports timely alerts to safeguard user health and safety.
- 隐私安全：无摄像隐私风险，在物理层和数据层完全严格遵守隐私法规，保护用户信息安全。
Privacy Security: No camera privacy risks, fully complies with privacy regulations at both physical and data layers to protect user information security.

1.3 应用场景 / Applications

- 商务办公：安装在办公室区域，实现人员在位监测，智能管理办公环境能耗与安全。
Business Office: Installed in office areas for occupancy monitoring, enabling intelligent management of office energy consumption and safety.
- 居家养老：安装在卧室或客厅，实时监测老人心率、呼吸和体动，异常情况及时警报。
Home Care: Installed in bedrooms or living rooms to monitor the elderly's heart rate, respiration and body movement in real time, with timely alerts for abnormalities.
- 商业场所：安装在酒店客房或公寓，实现入住检测与智能客房管理。
Commercial Places: Installed in hotel rooms or apartments for occupancy detection and smart room management.
- 学校宿舍：安装在宿舍区域，监测学生作息情况，保障住宿安全。
School Dormitories: Installed in dormitory areas to monitor students' daily routines and ensure accommodation safety.
- 月子中心：安装在母婴房间，监测产妇与婴儿体征，守护母婴健康。
Postpartum Care Centers: Installed in maternal rooms to monitor vital signs of mothers and infants, safeguarding their health.
- 康养机构：安装在养老院或康复中心，24 小时连续监测长者体征，提升照护效率。
Wellness Institutions: Installed in nursing homes or rehabilitation centers for 24/7 continuous vital signs monitoring, improving care efficiency.

1.4 外形结构 / Appearance Structure

1.4.1 外形尺寸 / Dimension


图 1.2: 产品外形尺寸图
Fig1.2: Product Outline Dimension Drawing

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***Note:** "Configurable?" means whether an interface is provided for users to set the parameter within the permissible range as required.

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The sensor should be placed between the mattress and the bed board. When a person is lying flat, the chest should be positioned directly above the sensor. A mattress thickness of 15~25 cm is recommended.

3.1.2 供电与电源要求 / Power Supply Requirements

- 使用 5V2A USB Type-A 接口的适配器进行供电。
Use a 5V2A USB Type-A adapter for power supply.
- 过大的电源纹波可能对雷达传感器造成干扰产生误报，建议供电电源纹波应低于 100mV。
Excessive power supply ripple may cause interference to the radar sensor and result in false alarms. It is recommended that the power supply ripple should be less than 100mV.

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Avoid objects that are swaying, moving, or significantly vibrating, as well as large metal surfaces within the radar detection range, including but not limited to: chandeliers, fans, air conditioners, plants, curtains, and mirrors with metal coatings.

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When installing the radar sensor, if there is an aluminum substrate or other metal plate on the back, it should be raised at least 5mm away from the metal surface and must not be in direct contact. The front of the radar sensor antenna should not be blocked by metal or covered by high-current cables. Avoid facing the drive power supply directly, and keep away from high-power components such as rectifier bridges, transformers, and switching transistors.

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The radar sensor has good penetration through plastic and wood materials, but cannot penetrate metal or metal-coated materials. If the product housing is made of special materials such as glass, ceramic, or carbon fiber, please refer to actual test results and contact EasyDetek technical personnel for applicability debugging if necessary.

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To achieve better presence detection and improve user experience, it is recommended to avoid movable objects within the detection coverage area, including but not limited to: fluttering curtains, swaying plants, rotating fans, etc.

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When the radar sensor is used together with wireless communication modules (NB, Bluetooth, WiFi, 2.4G modules), maintain sufficient spacing. It is recommended to keep a distance of at least 1m from high-power wireless communication devices such as routers and wireless hotspots during installation.

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The following statuses are reported in real time and support statistical analysis.

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Includes sleep onset time, wake-up time, sleep duration, sleep quality score, light sleep/deep sleep/REM ratio, and turning count.

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Includes presence/sleep status, body movement value, heart rate, respiration rate, heart rate status, and respiration rate status.

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The device is equipped with red and green LED indicators to indicate the current working status.

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3.3.2 查看与绑定设备 / View and Bind Devices

查看我的设备 / View My Devices

登录成功后可以看到账号的用户信息，点击导航栏【我的设备】即可看到该账号下绑定的所有设备。
After logging in, you can see the account user information. Click [My Devices] in the navigation bar to view all devices bound to the account.

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On the [My Devices] page, you can click the add button in the upper right corner to add (bind) a new device. Click any device card to enter the corresponding device details page, where you can view real-time monitoring data, historical data, sleep quality reports, etc. The above is the viewing process in the mini program; on the web portal, you can directly view all devices bound to the account after logging in.

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On the device details page, click the [1h / 6h / 12h / 24h] buttons at the top right of the chart to view data of different time periods; the sleep quality report can be viewed by clicking the [Today / This Week / This Month] buttons at the top right to view summary reports of different time periods.

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The settings interface allows you to view the current software version information of the device and configure parameters related to sleep and vital signs detection alerts.

To improve the accuracy of the sleep quality report, you need to set the time zone, the usual start time and end time of daily sleep. The daily sleep quality report is generated at the later of the set sleep end time and the actual sleep end time.

It is recommended to keep the detection distance range at the default value. If modification is required, please consult technical support personnel.

After confirming the configuration parameter changes, please click the Save Settings button.

3.3.5 体征实时状态与统计分析值 / Vital Signs Real-Time Status and Statistics

1. 体征实时状态值 / Vital Signs Real-Time Status:

- **在位/离位状态:** 当人进入检测区域后，即会显示在位，人离位一段时间后，会显示离位。
Presence/Absence Status: When a person enters the detection area, it will display "present"; after the person leaves for a while, it will display "absent".
- **清醒/睡眠状态:** 传感器检测到人睡眠后，分别会进入眼动、浅睡、深睡状态，并进行相应切换。
Awake/Sleep Status: After the sensor detects that the person is asleep, it enters REM, light sleep, and deep sleep states respectively and switches accordingly.
- **呼吸和心率值:** 当人进入检测区域后，不会立即显示呼吸和心率值，需等待被测对象处于静息状态后，才

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Respiration & Heart Rate Values: After a person enters the detection area, respiration and heart rate values are not displayed immediately; they are only displayed after the subject is in a resting state.

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Vital sign values can be output within as fast as 10 seconds after the person rests. If the person is not in a resting state, the sensor may not output valid vital sign values until the user is at rest.

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This product uses a USB interface (Type-C) for power supply and data communication. Connect the device to a 5V DC power adapter via a USB cable for power supply.

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This product consists of a radar monitoring device, a cloud service platform, and clients (App/Web). The communication system is divided into two layers: **device access** and **application interface**.

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The device establishes a persistent connection with the cloud via the MQTT protocol for data reporting, remote configuration, and firmware upgrade (OTA).

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- 数据上报 / Data Report: radar/{device_type}/{device_id}/data

样例 / Example: Topic:radar/edv28a/866965083771979/data, QoS: 0

```
{  
  "time_stamp": 1786333856,  
  "breath_state": 0,  
  "breath_confidence": 2,  
  "distance": 29,  
  "breath_rate": 15,
```

```
"sleep_state": 1,  
"heart_confidence": 1,  
"heart_rate": 52,  
"heart_state": 0,  
"movement": 10  
}
```

- 配置下发 / Config Set: radar/{device_type}/{device_id}/config/set

样例 / Example: Topic:radar/edv28a/866965083771979/config/set QoS: 0

```
{  
  "apnea_sec": 5,  
  "detect_range_max": 80,  
  "detect_range_min": 20,  
  "device_id": "866965083771979",  
  "device_name": "EDV28A_ 测试",  
  "hr_high": 120,  
  "hr_low": 50,  
  "interest_sleep_end_time": 360,  
  "interest_sleep_start_time": 0,  
  "no_vital_sec": 10,  
  "report_interval": 5,  
  "rr_high": 40,  
  "rr_low": 10,  
  "time_zone": 8  
}
```

- 配置确认 / Config Ack: radar/{device_type}/{device_id}/config/ack

样例 / Example: Topic:radar/edv28a/866965083771979/config/ack, QoS: 0

```
{  
  "result": "success",  
  "msg": "Config updated"  
}
```

- 固件请求 / OTA Request: radar/{device_type}/{device_id}/getOta

样例 / Example: Topic:radar/edv28a/866965083771979/getOta QoS: 0

```
{  
  "deviceType": "EDV28A",  
  "moduleType": "Air780EHV"  
}
```

- 固件响应 / OTA Response: radar/{device_type}/{device_id}/otaRes

样例 / Example: Topic:radar/edv28a/866965083771979/otaRes QoS: 0

```
{  
  "filename": "EDV28A-Air780EHV-v2.0.7-ota-0xF3620EF1.bin",  
  "ossUrl": "http://120.24.115.60:8080/firmware-files/EDV28A/Air780EHV/  
EDV28A-Air780EHV-v2.0.7-ota-0xF3620EF1.bin",  
  "status": "succse"  
}
```

4.2.2 应用接口协议 / Application Interface Protocol (REST/HTTPS)

客户端基于 HTTPS 的 RESTful API 与云端交互，提供用户管理、设备绑定、数据查询与配置管理。
Clients interact with the cloud via HTTPS-based RESTful APIs for user management, device binding, data query, and configuration management.

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Transport Encryption: Communication between device and cloud, as well as between client and cloud, is transmitted through encrypted channels.

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Authentication: Devices use username/password authentication; clients use Bearer Token authentication, with passwords encrypted via RSA for transmission.

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Data Isolation: Devices are bound to user accounts, and data is isolated by account.

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RF Installation Specifications: This device is equipped with a millimeter-wave radar module. Do not disassemble the device housing without authorization. During installation, avoid placing the sensor back tightly against large metal surfaces. There should be no fixed metal obstructions in front of the radar detection area to prevent affecting detection performance.

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Cable Usage Specifications: Do not pull or repeatedly bend the device's connecting cable to avoid internal damage.

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Environmental Requirements: The device must be operated and stored within the specified temperature range. Avoid prolonged exposure to humid or corrosive environments.

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Installation Layout Specifications: Keep the installation site away from electromagnetic interference sources such as water pumps, fans, high-power frequency converters, and wireless routers. Avoid persistent moving obstructions within the detection field of view.

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Usage Prohibitions: Do not modify the device hardware without authorization. Do not intentionally block the radar detection area. Do not use the device for information collection scenarios prohibited by laws and regulations.

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Disposal: Do not discard the device arbitrarily at the end of its service life. Please dispose of it in accordance with local electronic waste recycling regulations.

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Disclaimer for Improper Use: If the device is installed or used without following the safety compliance requirements in Section 5.1, resulting in device damage or abnormal detection function, the warranty shall be voided, and the company shall not be liable for any resulting losses.

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Disclaimer for Detection Limitations: Heart rate, respiration, presence/absence, and sleep status are identified through chest micro-movements. Reasonable deviations in monitoring data due to mattress conditions, blanket thickness, multiple persons in bed, or external vibration interference are considered normal.

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Disclaimer for Data Usage: All vital signs and sleep statistics output by the device are for daily health reference only. They shall not be used as a basis for professional medical diagnosis, and health conditions should not be determined solely based on the device monitoring results.

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Disclaimer for Network and Cloud Services: Data display on the Mini Program and Official Account platforms relies on mobile communication networks and cloud services. The company shall not be liable for data delays or loss caused by network failures or server maintenance.

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Disclaimer for Force Majeure: Device failures caused by external factors such as lightning strikes, abnormal power surges, strong electromagnetic interference, or natural disasters are not covered by the warranty.

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Disclaimer for Human-Induced Damage: Device damage caused by water ingress, external impact or extrusion, unauthorized disassembly or modification, or non-official repairs is not covered by the warranty.

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Version Iteration Note: The product hardware and software are continuously iterated and optimized. This manual is for reference only. Function adjustments and updates will not be separately notified.

5.3 售后服务 / After-Sales Service

图 5.1: 微信公众号
Fig5.1: WeChat Official Account


图 5.2: 微信小店
Fig5.2: WeChat Store